Medikoel

Me120

SI

APARAT ZA AEROSOLNO TERAPIJO
z ultrazvokom
Navodila za uporabo

GB

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Electromagnetic Compatibility Information

Table 1 - For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emissions		
The nebulizer Me120 is intended for use in the electromagnetic environment specified below. The customer or the user of the nebulizer Me120 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The nebulizer Me120 uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The nebulizer Me120 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Table 2 - For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity			
The nebulizer Me120 is intended for use in the electromagnetic environment specified below. The customer or the user of the nebulizer Me120 should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

Table 3 - For ME EQUIPMENT and ME SYSTEM that is not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity			
The nebulizer Me120 is intended for use in the electromagnetic environment specified below. The customer or the user of the Me120 should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2.5 GHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the nebulizer Me120, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Where the P is the maximum power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ¹ should be less than the compliance level in each frequency range. ² Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a- field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the nebulizer Me120 is used exceeds applicable RF compliance level above, the nebulizer Me120 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the nebulizer Me120. b- over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V] V/m.			

Table 4 - For ME EQUIPMENT and ME SYSTEM that is not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the nebulizer Me120	
The nebulizer Me120 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer of the user of the nebulizer Me120 can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the nebulizer Me120 as recommended below, according to the maximum output power of the communications equipment.	

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